

10k[®]

The World's Friendliest Pump

■ Easy to Load ■ Quick to Prime ■ High Performance



 **CONMED**[™]
L I N V A T E C

10k[®] Pump

The Simplest High-Performance Pump in Arthroscopy.

Ease of use and superior performance are the key characteristics that make the 10k Fluid System unique.

EASE OF USE

- Intuitive one-step/one-hand cassette loading
- Quick setup and priming
- Integrated pole mount works with any standard I.V. pole
- Compatible with video cart shelf or any flat surface

SUPERIOR PERFORMANCE

- Immediate response
- High flow
- Maintains distention under even the most demanding situations
- Robust autoclavable remote control

SPECIFICATIONS


- Maximum set pressure 150 mmHg
- Maximum flow rating 2000 ml/min
- Weight: 7 pounds
- Dimensions: 6.5" (16.5 cm) tall, 8" (20.3 cm) deep and 7.75" (19.7 cm) wide
- Voltage: 90 – 240 v



10k Console



C7115 Autoclavable Remote Control



10k100 Arthroscopy Inflow Tube Set



10k150 Arthroscopy Inflow/Outflow Tube Set



10K600 Eco-Flow™ Day-use Tubing Set



10K605 Eco-Flow™ Single-use Patient Tubing

ORDER INFORMATION

Cat. No.	Description
10k	Console
10k100	Arthroscopy Inflow Tube Set
10k100D	Arthroscopy Inflow Tube with Colder Set
10k150	Arthroscopy Inflow/Outflow Tube Set
10k505	Laparoscopic Irrigation/Suction Tube Set
10k600	Eco-Flow™ Day-Use Pump Tubing <i>(Must be used with 10k605)</i>
10k605	Eco-Flow™ Single-use Patient Tubing <i>(Must be used with 10k600)</i>
10k009	(Optional) Console Mounting Kit
C7115	Autoclavable Remote Control

ConMed Linvatec Fluid Management System

11311 Concept Boulevard
Largo, Florida 33773-4908

Local 727-392-6464
Toll Free 800-237-0169

Website www.linvatec.com
Email customer_service@linvatec.com