



High Definition Scopes and QuickLatch®

- 4.0mm HD Scope Systems
- Quicklatch Sheaths (4.0mm scopes)
- True autoclavability
- Eyepiece and Cartridge available



CLEARLY TRUE HIGH DEFINITION™

True HD



High Definition Scopes

ConMed Linvatec brings another first to market with the HD arthroscope. Continuing our drive of optical optimization, the enhanced lens design ensures maximized transmission of high contrast light from the arthroscope into the True HD camera head. ConMed Linvatec can now offer the highest quality optics throughout the whole Imaging chain for a True HD image in ENT and Orthopedics. Our HD scope range is compatible with your existing Quicklatch® sheath range ensuring you don't have to settle for low flow to operate in high definition.


- Improved light transmission
- High contrast to the picture edge
- Improved depth of field
- High contrast to the back of the joint

	EYECUP	CARTRIDGE
0 Degree	HD4000	HD4001
30 Degree	HD4300	HD4301
70 Degree	HD4700	HD4701
30 Degree Rev	HD4302*	

*For use with QLR6131 or QLR6132 sheath sets only.



Cartridge



Sheath Systems

QuickLatch® 4.0mm Sheaths

QuickLatch design increases durability and provides excellent distal tip protection. Tangent position of scope and proprietary stopcock design provide excellent flow rates in both the 5.5mm and 5.9mm sizes

	SINGLE	DUAL
5.5mm Arthroscopy	QL5531	QL5532
5.9mm High Flow	QL6131	QL6132
5.9mm Two-Connect	QL6231	QL6232

See Back Panel for Trocar and Accessory Items



Bridge Cannulas



Multiple sizes for multiple requirements. Cannulated obturators are compatible with 3.175mm switching sticks and portal rods.

	4.6mm	5.8mm	7.6mm
	Inner Diameter		
Standard	OC4601	OC5801	OC7601
Standard 30/70°	OC4606	OC5806	OC7606
Trocar	OC4603	OC5803	OC7603
Conical Obturator	OC4604	OC5804	OC7604
Cannulated Obturator	OC4605	OC5805	OC7605

QuickLatch 4.0mm Bridge System

The QuickLatch Bridge for 4mm scopes is available in both single and double stopcock versions. Two-Connect adapter allows for "bridging" with ConMed Linvatec's pressure sensing pumps (4.6mm Cannulas are not compatible).



Single Stopcock	QL6111
Dual Stopcock	QL6112
Two-Connect	QL6210*

*With QL6230, T5101 or T5102